

Allfit™

GENERAL PRODUCT INFORMATION

Allfit are low temperature thermoplastic pellets that can be kneaded and formed to any shape. This product is used in physical rehabilitation to make individual web spacers, add-ons to orthoses, to fabricate customized handles, utensils, pencils and all kinds of tools. The pellets can also be used for desensitization exercises.

PRODUCT RANGE

Article Number	Description
703011	Allfit™

PRECAUTIONS BEFORE USE

- High temperatures up to a maximum of 120°C (248°F) do not damage Allfit pellets, but are not user-friendly.

WORKING TECHNIQUE

1. Pour Allfit in a Suspan, a water heater or a bowl of hot water at 65°C (149°F).
2. The pellets will stick together and become transparent to form a viscous mass.
3. Knead the Allfit pellets homogeneously.
4. Form it around an object into the desired shape. Allfit sticks to most orthotic fabrication materials that are not coated and to porous surfaces.
5. Cooling can be enhanced by means of cold water or cold spray.
6. Allfit is completely or partially remouldable by reheating the moulded Allfit in hot water of 65°C (149°F) it is possible to correct and reshape the Allfit into the desired form or to remove it from the chosen object and apply onto another object.

Beware: Temperatures of 65°C (149°F) or more can also be reached in the patient's daily life. Think of a closed car in the summer, the surface of a hot radiator, a sauna, or the proximity of an open fireplace.

WASTE MANAGEMENT

Allfit can be disposed of with normal household waste without harming the environment.

STORAGE

- It is advised to store the product in a dark, cool, dry place in the original packaging at temperatures of min. 10°C (50°F) and max. 30°C (86°F).
- Once removed from the packaging, the leftovers should be repackaged to avoid biodegradation.
- Do not lose the labels. They are crucial for product traceability.

GENERAL SAFETY ADVICE

- ***This product is not suitable for internal use. It may not be used on open wounds or in the mouth.***
- ***Never use an open flame to activate the material.***
- ***orthoses and rehabilitation aids made out of Orfit thermoplastics can only be made by qualified health professionals.***
- ***In the event of a serious incident that is related to the use of this product, you must report this to the distributor in your country. Please consult our website to find your distributor.***

ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets, and regulatory information, please visit our website www.orfit.com.

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.
This product is a registered trademark from ORFIT INDUSTRIES BV.



ORFIT INDUSTRIES BV

Vosveld 9A | B-2110 Wijnegem | Belgium

T (+32) (0)3 326 20 26

welcome@orfit.com

www.orfit.com



Ref. Nr. 31150

Version 32

Last updated 1/12/2024

Revision date 1/12/2026