

## LUXOFOAM

### GENERAL PRODUCT INFORMATION

LUXOFOAM is a self-adhesive, closed-cell soft PVC foam.

It is intended for use as cushioning and padding material in applications such as the fabrication and adjustment of custom-made support, positioning, and protective solutions.

### PRODUCT RANGE

LUXOFOAM is available in sheets of 50 cm (width) x 97 cm (length) x 3 mm (thickness).

| Art. no. | Perforation type |
|----------|------------------|
| 35110    | non-perforated   |
| 35112    | perforated       |

### PRECAUTIONS BEFORE USE

- Ensure the workplace is well ventilated.
- Use appropriate tools and ensure safe handling during preparation and assembly.
- Ensure a comfortable and stable working position during use.
- Cut the required pattern before removing the protective paper from the adhesive backing.

### WORKING TECHNIQUE

- Remove the protective paper from the adhesive backing and apply LUXOFOAM to the desired surface.
- LUXOFOAM may be applied before activating the fabrication material or after completion of the fabricated item. When applied afterwards, consider the additional thickness added by the foam.
- LUXOFOAM can be repositioned or removed easily immediately after application. Removal may become more difficult over time unless the applied surface is reheated according to the properties of the underlying material.

### MAINTENANCE

- LUXOFOAM can be cleaned with lukewarm water and liquid soap. Rinse thoroughly and dry completely.

**⚠ If you are not sure about the cleaning agent, do not use it.**

**⚠ Do not use solvents or abrasive detergents.**

- LUXOFOAM is not suitable for sterilization in an autoclave.

### WASTE MANAGEMENT

LUXOFOAM can be disposed of with normal household waste after use.

### USER INFORMATION

The qualified professional responsible for the preparation, application, or integration of the product into the final application shall provide the user with appropriate information regarding its use, maintenance, and limitations. The user should follow the instructions provided for the correct use of the final application. Contact the qualified professional if adjustments, modifications, or further guidance are required.

### STORAGE

- It is advised to store the product in a dark, cool, dry place in the original packaging at the temperature of min. 10°C (50°F) and max. 30°C (86°F).
- Do not lose the labels, They are crucial for product traceability.

### GENERAL SAFETY ADVICE

- For external use only. Do not use it on open wounds or in the mouth.
- LUXOFOAM may only be modified or incorporated by qualified professionals.

### ADDITIONAL INFORMATION

For additional information, such as distributor contact information, product brochures, safety data sheets, and regulatory information, please visit our website at [www.orfit.com](http://www.orfit.com)

#### Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.

#### ORFIT INDUSTRIES

Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

[www.orfit.com](http://www.orfit.com)



Ref. No. 31300

VERSION 6

LAST UPDATE: 13/07/2026

REVISION DATE: 13/07/2028