

ORFILASTIC

GENERAL PRODUCT INFORMATION

ORFILASTIC is an elastic terrycloth lining material with a self-adhesive open-cell foam backing. It is intended for use as a lining and cushioning component in custom-made applications where improved air permeability and surface comfort are desired.

PRODUCT RANGE

ORFILASTIC is available in sheets of 50 cm (width) x 100 cm (length) x 3 mm (thickness).

Art. No.	Colour
35111	Beige
35117	Black

PRECAUTIONS BEFORE USE

- Ensure the workplace is well ventilated.
- Use appropriate tools and ensure safe handling during preparation and assembly.
- Ensure a comfortable and stable working position during use.

WORKING TECHNIQUE

- Cut the required pattern before removing the release paper from the adhesive backing. The cutting direction has no influence on product performance.
- Apply ORFILASTIC to the required surface after the final shape of the application has been obtained.
- Remove the release paper gradually while applying ORFILASTIC to the surface with slight pressure. Avoid creating folds or wrinkles during application.
- Ensure that ORFILASTIC does not come into contact with itself (adhesive side to adhesive side), as separation is not possible afterwards.

MAINTENANCE

ORFILASTIC is NOT waterproof and has to be replaced when necessary.

WASTE MANAGEMENT

ORFILASTIC can be disposed of with normal household waste after use.

USER INFORMATION

The qualified professional responsible for the preparation, application, or integration of the product into the final application shall provide the user with appropriate information regarding its use, maintenance, and limitations. The user should follow the instructions provided for the correct use of the final application. Contact the qualified professional if adjustments, modifications, or further guidance are required.

STORAGE

- It is advised to store the product in a dark, cool, dry place in the original packaging at the temperature of min. 10°C (50°F) and max. 30°C (86°F).
- Do not lose the labels, They are crucial for product traceability.

GENERAL SAFETY ADVICE

- For external use only. Do not use it on open wounds or in the mouth.
- ORFILASTIC may only be modified or incorporated by qualified professionals.

ADDITIONAL INFORMATION

For additional information, such as distributor contact information, product brochures, safety data sheets, and regulatory information, please visit our website at www.orfit.com

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries.

ORFIT INDUSTRIES

Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. No. 31301

VERSION 7

LAST UPDATE: 13/07/2026

REVISION DATE: 13/07/2028