







AIO Solution 3.0 - Belly & Pelvic Board Components

Individual patient components list

Patient ID:

Date:

	YES Art. No.:	NO		38126	38127
AIO 3.0 Base Plate			38128		
	38101 - Adult			38204	
	38117 - Paediatric			38210	
AIO 3.0 - 0° cushion			Extension frame for knee & foot support		
	38107			38206	
AIO 3.0 - Base cushion			Indexable Foot Support		
	38108 - Soft Chest Cushion			38226	
AIO 3.0 - Insert Cushion			AIO 3.0 - Indexed Prone Head Support		
	38109 - Full Insert Cushion			Model 1	Model 2
				Model 3	Model 4
				Model 5	Model 6
AIO 3.0 - Insert Cushion			High precision head supports	Low density	Regular density
	38204			20 mm	40 mm
	38210			9°	13,5°
Extension frame for knee & foot support				18°	
			High precision blocks & wedges	Low density (CFL)	Low density (GFL)
				High density	

All rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries ©ORFIT Industries

for more information please visit our website at www.orfit.com

Art. No. 22156E - V2 - 27/08/2024





AIO Solution 3.0 - Breast- & Lungboard Components

Individual patient components list

Patient ID:

Date:



Positioning of elements									
Position of the patient			Prone			Supine			
AIO base plate position (position of laser line on millimeter scale):									
Set-up without base plate: 2-pin bar positions for AIO Base Cushion			Upper			Lower			
Position of the knee & leg immobilisation cushions on the extension frame			Frame 38204 (1-19)			Frame 38210 (1-27)			
Side of the knee & leg immobilisation cushions on the extension frame			A			B			
Position of the knee & leg immobilisation cushions on the extension frame			Frame 38204 (1-19)			Frame 38210 (1-27)			
Position of the indexable foot support on the extension frame			Frame 38204 (1-19)			Frame 38210 (1-27)			
Indexable foot support adjustments			Tilt position			Foot cup rotation			
			1	2	3	A	B	C	D
Position of the knee & leg immobilisation cushions on a 2-pin bar			1-5			2-pin bar position			
Position of head support			1			2		3	
Leg separator position			Cranial			Caudal			
Thermoplastic masks used:			4-points pelvic mask prone/supine						
			6-points pelvic mask prone/supine						
			6-points pelvic mask supine						
Others/remarks									

All Rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries ©ORFIT Industries