



Patient ID:

Date:


	YES Art. No.:	NO
AIO 3.0 Base Plate		


	38101 - Adult
	38117 - Paediatric
AIO 3.0 - 0° cushion	

	38110
AIO 3.0 - Cranial Cushion	


	38111
AIO 3.0 - Caudal Cushion	

	38126	38127
	38128	
Knee & leg immobilisation cushions		

	38204
	38210
Extension frame for knee & foot support	


	38112
AIO 3.0 - Cranial Elevation Cushion	

	38113
AIO 3.0 - Caudal Elevation Cushion	

	38114
AIO 3.0 - Comfort Wedge for Contralateral Breast	

	38115
AIO 3.0 - Bridge Plate for Contralateral Breast	

	38219
AIO 3.0 - Indexation Cushion for use without Base Plate	

	38226
AIO 3.0 - Indexed Prone Head Support	29051
	Comfort Head Rest for Prone Position

All Rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries ©ORFIT Industries



AIO Solution 3.0 - Prone Breast Board Components

Individual patient components list

Patient ID:

Date:

	35703 (low density)	35706 (low density)
	35703/MR (low density)	35706/MR (low density)
	35754/1 (high density)	35754/2 (high density)
Head Support Elevation Blocks (20 or 40 mm)		

	Longitudinal position					
	1	2	3	4	5	6
	Rotational position					
	A		B		C	

AIO 3.0 - Hand Grip for Prone Breast Board

Positioning of elements			
AIO base plate position (position of laser line on millimeter scale):			
Set-up without base plate: 2-pin bar position for AIO 3.0 Cranial Cushion			
Set-up without base plate: 2-pin bar position for AIO 3.0 Indexation Cushion for use without Base Plate			
Position of the left hand on the AIO 3.0 Cranial Cushion	1	2	3
	4	5	6
Position of the right hand on the AIO 3.0 Cranial Cushion	1	2	3
	4	5	6
Position of the Indexed Prone Head Support	Position 1		Position 2
Treated Breast	Left		Right
Direction of the Comfort Wedge for Contralateral Breast	Highest point medial		Highest point lateral
Side of the knee & leg immobilisation cushion on the extension frame	A		B
Position of the knee & leg immobilisation cushions on the extension frame	Frame 38204 (1-19)		Frame 38210 (1-27)
Position of the knee & leg immobilisation cushions on a 2-pin bar	1-5		2-pin bar position
Thermoplastic mask used:	Model:		
Others/remarks			

All Rights reserved. No part of this publication may be produced, stored in a retrieval system or transmitted, in any way, form or by any means, electronic, mechanical, photocopying or otherwise, without prior written permission of Orfit Industries © ORFIT Industries