

Orfitrans is an advanced range of high temperature thermoplastics providing solutions for a wide variety of prosthetic and orthotic applications.

The Orfitrans product family is available in various configurations to accommodate different size needs, as well as patient and clinical preferences.

Orfitrans PETG

Orfitrans Stiff

Orfitrans Excel

Orfitrans Extra Soft

Orfitrans Extra Soft Silicone



The reliable choice

IN PROSTHETIC AND ORTHOTIC THERMOPLASTIC MATERIALS



ORFITRANS®

Thermoplastics for test sockets, inner sockets and orthotic devices



The reliable choice

IN PROSTHETIC AND ORTHOTIC THERMOPLASTIC MATERIALS

Orfit Academy

Orthotic & Prosthetic Institute

To learn more about Orfitrans products, visit the Orfit Academy, Orthotic and Prosthetic Institute, where we have an overview course relating to all Orfitrans orthotic and prosthetic fabrication materials.



Go to www.orfit.com for more product information, such as article numbers, instructions for use, and product education.

Visit
WWW.ORFIT.COM
for more information.

ORFIT INDUSTRIES
Vosveld 9A
B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcomel@orfit.com

ORFIT INDUSTRIES AMERICA
810 Ford Drive
Norfolk, VA 23523 | USA
T 516-935-8500
welcomel@orfit.com

08/07/2024
Art. N° 53002

The Orfitrans Product Line consists of materials for the fabrication of test sockets, inner sockets and orthotic devices.

Use this overview to select the right product for your next O&P device.

Product Specifications

Thickness		Sheet Size		Art. No. Orfitrans PETG		Art. No. Orfitrans Stiff		Art. No. Orfitrans Excel		Art. No. Orfitrans Extra Soft		Art. No. Orfitrans Extra Soft Silicone	
mm	Inch	mm	Inch	Transparent	Transparent	Transparent	Transparent	Semi-transparent	Black	Semi-transparent	White		
3	1/8	1225 x 1225	48 x 48	x	x	x	x	37SS42/1225122530	37SS42/1225122530Z	x	x		
5	13/64	1225 x 1225	48 x 48	x	x	x	x	37SS42/1225122550	37SS42/1225122550Z	x	x		
6	15/64	400 x 400	15 ¾ x 15 ¾	x	x	x	x	x	x	x	x		
8	5/16	406 x 406	16 x 16	x	37ST/4480	x	37ST/4480	37SS42/4480	37SS42/4480Z	x	x		
9	23/64	400 x 400	15 ¾ x 15 ¾	x	x	x	x	x	x	37SS62/4490	37SIL/4490		
10	3/8	406 x 406	16 x 16	37PETG/44100	37ST/44100	x	37ST/44100	37SS42/44100	37SS42/44100Z	x	x		
12	1/2	406 x 406	16 x 16	37PETG/44120	37ST/44120	x	37ST/44120	37SS42/44120	37SS42/44120Z	37SS62/44120	37SIL/44120		
15	19/32	406 x 406	16 x 16	37PETG/44150	37ST/44150	x	37ST/44150	37SS42/44150	37SS42/44150Z	x	37SIL/44150		
18	7/10	400 x 400	15 ¾ x 15 ¾	x	37ST/44180	x	37ST/44180	x	x	x	x		
20	3/4	406 x 406	16 x 16	37PETG/44200	x	x	x	x	x	x	x		



Orfitrans PETG

The clear appearance lets the prosthetist/orthotist verify the conformity of the device prior to final fabrication and after setting. The product has superior impact resistance and is exceptionally durable in handling, drilling, sanding, buffing, cutting and re-forming making its clinical application and adjustment effortless for clinicians.

Recommended application

Prostheses or orthoses for short-term, supervised, in-office test fittings



Orfitrans Stiff

Orfitrans Stiff is the Gold Standard for check sockets. The material is virtually indestructible, yet easy to drill into or modify. The socket will remain transparent with a slight blue tint that allows the limb to remain visible.

Recommended application

Prosthetic test sockets for long-term test applications.



Orfitrans Excel

8% lighter than competitive EVA-type products, resulting in increased patient comfort. The material's self-bonding capabilities make it capable of all fabrication techniques for both orthoses and prostheses.

Recommended application

Inner prosthetic sockets, hand and wrist orthoses, and as a lining material in ankle orthoses.



Orfitrans Extra Soft

Highly flexible, very stable product that is extra soft and smooth for added patient comfort.

Recommended application

Inner prosthetic sockets for very active amputees.



Orfitrans Extra Soft Silicone

Extra soft, smooth material with the most silky surface of all Orfitrans products. Thanks to the silicone, this product prevents the need for powder on the inside and it eliminates all the difficulties and inconveniences often experienced in donning and doffing other socket materials.

Recommended application

Inner prosthetic sockets for very active amputees.


The Orfitrans product family

Test Sockets

Orfitrans PETG for short-term, supervised test fittings




Orfitrans Stiff for long-term prosthetic test sockets



Inner Sockets

Orfitrans Excel



Orfitrans Extra Soft



Orfitrans Extra Soft Silicone



Orthotics

Orfitrans Excel for hand, wrist, and ankle lining orthoses

